



DENTAL-EVAC Kit

FSP6000

Kit contains:

- ◆ Reusable DENTAL-EVAC pump
- ◆ Two canisters with FSP filter
- ◆ Wide bore soft yankauer
- ◆ Saliva Ejector
- ◆ Carry Bag



DENTAL-EVAC Replacement

FSP6000R

Replacement contains:

- ◆ Two canisters with FSP filter
- ◆ Wide bore soft yankauer
- ◆ Saliva Ejector
- ◆ Sterile 14fr, 10fr and 8fr catheters

Sold in quantities of 10

DENTAL-EVAC®

Hand-Held Suction For **Dentists**

Easy to use as 1-2-3!



Select canister and suction tip



Snap canister into reusable pump



Introduce suction tip into patient's airway
squeeze handle at desired rate to suction

DENTAL-EVAC®

Hand-Held Suction For **Dentists**



- ✓ Portable, hand powered suction
- ✓ Ideal for emergency suction back-up
- ✓ Wall suction level performance in a portable unit >600 mmHg
- ✓ Always ready to go – no set up required, no batteries to recharge or replace
- ✓ Powerful enough to suction in all situations
- ✓ Patented 0.22 equivalent micron filter ensures maximum infection control
- ✓ Convenient one-handed operation
- ✓ Each unit is 100% tested for reliability



Sedation Resource

Keeping Sedation Simple

800-753-6376

P/N Rev. A



Sedation Resource

Keeping Sedation Simple

800-753-6376

DENTAL-EVAC®

Hand-Held Suction With FSP®

FSP® (FULL STOP PROTECTION) FILTER

- ✓ Blocks all known fluids and pathogens for total protection
- ✓ Protects DENTAL-EVAC pump for safe reuse

MEETS OR EXCEEDS

- ✓ OSHA 29 CFR1910.1030 Occupational Exposure to Bloodborne Pathogens
- ✓ CDC guidelines for SARS



DENTAL-EVAC®

Hand-Held Suction For Dentists

IDEAL FOR

- ✓ Emergency Back-Up Suction
- ✓ Surgery
- ✓ Anesthesia
- ✓ Retrieving a Lost Crown
- ✓ Lobby Emergencies



Technical Specifications

FSP® (Full Stop Protection) Filter

- Equivalent to 0.22 micron filter
 - Laboratory challenge tested with *Brevundimonas diminuta*
 - 100% Blockage of fluids, position insensitive
 - Blocks against the aerosolization of pathogens and viruses
- *Laboratory report available on request



DENTAL-EVAC Suction Pump

- Weight 143 grams (5oz)
- Size 114 mm x 180 mm x 55 mm (4.49" x 7.10" x 2.15")
- Volume per stroke 30 ml
- Peak pump vacuum >600 mmHg (>24" Hg)
- Peak air flow 20 L (50 ms)
- Internal Seals – flourosilicone (pump can be sterilized via ETO or Sterad)
- Reusable



Quality Assurance Level

- 100% testing of pump performance
 - 100% testing hermetic seal every canister and cap*
 - All DENTAL-EVACs built under Class 100,000 controlled environment conditions for known low bioburden
- *Quality Conformity Certificate available on all products



Canister Size

- 300 ml
- L152.4 mm (6")



Catheters Available

- Wide bore soft yankauer ID 9.53 mm (.375") L 272 mm (10.7")
- Saliva ejector tip ID 4.7 mm (.184") L 152.4 mm (6")
- Saliva ejector tubing ID 5.7 mm (.223") L 635 mm (25")



Protects Healthcare Providers from Pathogens and Viruses
Including Flu, Hepatitis, TB, HIV/AIDS, SARS and Avian H5N1 Bird Influenza

CALL 800-753-6376 FOR MORE INFORMATION