



iM3
Vital Signs Monitor

A Slim Monitoring Experience

iM3 Vital Signs Monitor



About Edan

Edan is a healthcare company dedicated to improving the human condition around the world by delivering value-driven, innovative and high-quality medical products and services. For over 20 years, Edan has been pioneering a comprehensive line of medical solutions that address a broad range of healthcare practices including:

- Diagnostic ECG
- Patient Monitoring
- OB/GYN
- Ultrasound Imaging
- Point-of-Care Testing
- In-Vitro Diagnostics
- Veterinary

Healthcare professionals around the world depend on Edan's breakthrough medical technologies and outstanding customer support.

Authorized Distributor:
Sedation Resource, Inc
800.753.6376

info@sedationresource.com



A world of potential

© 2017 Edan Instruments, Inc. All rights reserved. Features and specifications are subject to change without prior notice. No reproduction, copy or transmission may be made without written permission. Not all products or features are available in all countries, contact Edan for local availability.



ENG-PM-iM3-V1.0-20170928



A world of potential

Smarter, Larger, yet Slimmer

Combining a stunning ultra-slim design, a large capacitive touch screen, together with a rich set of leading functions, iM3 takes the lead and presents itself as the beginning of a new generation of vital signs monitors. Bringing in full options of SpO₂, NIBP and Quick Temp/Infrared Ear Temp and differentiated work modes, iM3 serves for both long time monitoring and spot check applications.



8" Touch Screen



360° Alarm Light




Cable-receiving Design




<2.5 kg

Ultra-thin & Light Weight



Built-in Barcode Scanner



TEMP Options

Standard Parameters

EDAN SpO₂, EDAN NIBP, PR

Optional Parameters

Nellcor™ SpO₂, SunTech® NIBP, EDAN Quick TEMP, Infrared Ear TEMP, Covidien Quick TEMP

Redefinition of Adaptive Work Modes

Round Mode

Special work mode designed for daily ward rounding application. Worklist could be easily input through USB disk or downloaded from EMR system. Moreover, the color-coded patient list comes with square windows interface design, you can tell patients' status and information in a quick glance.



Patient management



Patient status with color code

Monitoring Mode

For continuous patient monitoring with real-time data, real-time alarms, and history trend.

Spot Mode

For spot-check applications, such as out-patient, emergency, etc. Real-time vitals are acquired and displayed, no alarm will be seen or heard.

TEMP Module Options

TEMP NOT INCLUDED



Infrared Ear TEMP



EDAN Quick TEMP



Covidien Quick TEMP



Rolling Stand

Non-volatile built-in Memory

Spot & Round Mode

16,000,000

Measurements

Monitoring Mode (for each patient file)

200

Alarm

1200

NIBP Measurements

120 h

Trend