

Parker Reinforced Laryngeal Mask Device™

DISPOSABLE SILICONE Laryngeal Mask

Parker's reinforced disposable Silicone Mask Device is gentle to the airway.

Spiral-wound reinforcing wire in the Parker flexible silicone tube reduces the risk of kinking.

The smooth contouring of our reinforced laryngeal mask device enables easy insertion and is designed to produce an effective seal.

The cuff's matt finish is designed to reduce the risk of the mask slipping out of place.

Our disposable mask eliminates costly reprocessing and the risk of cross-contamination.

The unique color-coded and printed pilot balloon allows for fast and correct size identification.

The clear 15mm connector allows potential blockages at that site to be seen.

Parker's Laryngeal Mask Devices are conveniently packaged: five units per shelf box.

Cost effective
Latex-free
DEHP-free

When Softness Counts Count on Parker

Order

H-PLMDR-10	Size 1.0
H-PLMDR-15	Size 1.5
H-PLMDR-20	Size 2.0
H-PLMDR-25	Size 2.5
H-PLMDR-30	Size 3.0
H-PLMDR-40	Size 4.0
H-PLMDR-50	Size 5.0

All product is packaged 5/Box, 50/Case



To Life!™



Parker Medical, Inc.
9457 S. University Blvd.
PMB 331
Highlands Ranch, CO 80126
Phone: 303-799-1990
FAX: 303-799-1996